

Canon



Infinix-i

System configurations to
meet your clinical needs

Infinix-i Sky

R

C

Unprecedented flexibility and integration

The Infinix-i Sky is designed to move around you and your patient so procedures can be performed more comfortably, safely and efficiently.


- Achieve optimal angulations with unprecedented head-to-toe and fingertip-to-fingertip coverage
- Continuously maintains a heads-up display during compound angles with synchronized rotating collimator and FPD
- Access halo ensures unobstructed head-end work space to improve patient access for staff and ancillary equipment

Available detector sizes

20 x 20

30 x 30

30 x 40





270° C-ARM
ROTATION

Canon

CANON MEDICAL SYSTEMS EUROPE B.V.

129626, Москва, проспект Мира,
дом 102, корпус 1, этаж 6, к. 6

Телефон: 8 (800) 555-73-87

Email: info@medeq.ru

Web: www.medeq.ru

Infinix-i

Rite Edition

Made For life